



**ACCUMAX**  
SAFETY

# Plum Body Neutraliser

**Product Code:** AM-H4735

The Plum Body Neutraliser is a mobile body shower that is an effective decontamination solution for skin, surfaces and clothing.

For fast and effective decontamination of bodies and surfaces in contact with corrosive and toxic hazardous substances such as strong acids and bases, oxidants, formaldehyde, phenols, organic solvents, hydrofluoric acid and fluorides.

## Product Benefits

- For first aid after contact with dangerous substances directly at the scene of the accident
- Very rapid effect
- Top quality mobile compression sprayer with 5L Plum Body Neutral solution
- Easy handling: directly applicable after 10 pumping strokes
- With handle and stainless steel lance
- For decontamination of several persons or areas up to 50 m<sup>2</sup>
- Container made of lightweight, chemical-resistant special plastic



### Avoid eye contact!

In case of eye contact, rinse immediately with eye wash solution and seek medical advice. If possible show the label.



**ACCUMAX**  
GLOBAL

1300 222 862 | [sales@accumax.com.au](mailto:sales@accumax.com.au)  
58 Catalano Circuit, Canning Vale WA 6155  
[accumaxglobal.com.au](http://accumaxglobal.com.au)

## 5L Plum Body Neutral Solution

- Easy handling: directly applicable after 10 pumping strokes
- With handle and stainless steel lance
- For decontamination of several persons or areas up to 50 m<sup>2</sup>
- Spraying time: 3 min
- Shelf life: 3 years from date of manufacture
- Container made of lightweight, chemical-resistant special plastic
- Dimensions:  
Height: 570 mm, with lance 800 mm  
Width: 220 mm  
Diameter: 180 mm



### Plum Body Neutral Owes its Effectiveness to these Substances

- Amino acid buffer for neutralizing acids and alkalis such as sulphuric acid, caustic soda, nitric acid or hydrochloric acid
- Amino acid glycine for binding formaldehyde and phenols
- Calcium chloride for precipitation of fluoride as insoluble calcium fluoride, thereby also inactivation of hydrofluoric acid
- Additionally, polyethylene glycol and glycerol as solubilizers for organic chemicals

### Permits/Approvals of the Decontamination Solution

- Certification as a class I medical device, according to the requirements of the EU Regulation EU 2017/745
- Officially registered by the Federal Office for Drugs and Medical Devices under the number DE/CA33/ILM/2020/1. The registration is valid for the entire EU and may be marked with the CE mark accordingly
- Conforms to DIN 14555-12 for hazardous goods equipment trolley

